

One sample. Multiple results. Excellent sensitivity and specificity.

Aptima® Performance Data on the Panther® System

Specimen Type	Chlamydia ¹		Gonorrhea ¹		Trichomonas ²	
	Sensitivity	Specificity	Sensitivity	Specificity	Sensitivity	Specificity
Vaginal Swab (Clinician-Collected)	97.2%	98.5%	97.7%	99.6%	100%	98.20%
Vaginal Swab (Patient-Collected)	97.2%	95.5%	97.7%	99.6%	*	*
Female Endocervical Swab	97.2%	99.3%	100.0%	99.8%	100%	98.10%
Male Urethral Swab	100.0%	99.1%	100.0%	100.0%	*	*
Urine (Male)	95.2%	99.8%	98.7%	99.7%	*	*
ThinPrep® Liquid-based Pap	98.8%	100.0%	100.0%	100.0%	100.0%	98.60%

*These sample types are not FDA-cleared for Trichomonas testing.

1. Gen-Probe, 2013 APTIMA COMBO 2 Assay (package insert) #502446, Rev. 001, San Diego, CA, 2013
 2. APTIMA Trichomonas vaginalis assay (package insert) #503684, Rev. A, San Diego, CA

TRN-00117-001 Rev. 001 HOLOGIC, Aptima, Panther, ThinPrep, and Tigris are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries.

HOLOGIC